

## ISOTOPIC SEAWATER ANALYSIS USING CAVITY RING-DOWN SPECTROSCOPY: METHOD FOR REDUCING THE EFFECT OF SALTS ON THE MEASUREMENT


David Kim-Hak <sup>(1)</sup>, Kuan Huang <sup>(1)</sup>, Renato Winkler <sup>(1)</sup>, Claire Normandeau <sup>(2)</sup>, Elizabeth Kerrigan <sup>(2)</sup>, D.W.R Wallace <sup>(2)</sup>

<sup>(1)</sup> Picarro, Inc., 3105 Patrick Henry Drive, Santa Clara, California

<sup>(2)</sup> Department of Oceanography, Dalhousie University, Halifax, Nova Scotia, Canada

### Introduction


- Stable isotopes of water,  $\delta^{18}\text{O}$  and  $\delta^2\text{H}$ , are unique tracers for studying hydrological and associated processes.
- CRDS analyzer produces high-precision measurements of water isotopes in the gas phase.
- High total dissolved solids (TDS) waters can foul the vaporizer affecting the memory performance and requiring more frequent cleaning and greater downtime<sup>[1]</sup>.
- The Salt Liner is an accessory that protects the analyzer system from fouling while maintaining optimum performance



### CRDS Analyzer: L2130-i and L2140-i.

#### Cavity Ring-Down Spectroscopy (CRDS) Technology:

- Laser-based technique to measure the stable isotopes of water,  $\delta^{18}\text{O}$  and  $\delta^2\text{H}$ .
- Time-based measurement rather than absolute absorbance.
- Robust, compact and field-deployable.
- Highest precision and lowest drift.



### Salt Liner (SL)

- Stainless Steel mesh insert for the High-Precision Vaporizer(A0211)
- Catches salt precipitates as the injected High-TDS water sample is evaporated in the Vaporizer chamber
- Cleanable and reusable



### Test Method


- Precision test: 160 injections of a sample
- Memory test: 76 injections per sample. Two samples of different isotopic composition used.
- Sample to Sample Measurements: alternating 10 injections of two samples

#### - Sample Solutions:

- DI water
- Kona Deep Bottle Water
- Instant Ocean mixed with DI water and Kona Deep (4% Salinity)

### Results


#### 1) Drift, Precision, and Memory Test of DI water: with SL vs without SL



	Specifications		Measured w/o SL		Measured with SL	
	d18O	dD	d18O	dD	d18O	dD
Precision (permils)	>0.025	>0.1	0.01	0.05	0.01	0.05
Drift (permils)	>0.2	>0.8	0.06	0.33	0.06	0.2

- It takes 4 injections to get rid of 99% of the memory in  $\delta^{18}\text{O}$  and 98% memory in  $\delta\text{D}$
- Measurement performance of DI water is unaffected by the SL insert as Precision, Drift and Memory performance remain unchanged.

#### 2) Drift, Precision, and Memory Test of 4% Seawater Solutions with SL




	Specifications (with DI)		Measured with SL	
	d18O	dD	d18O	dD
Precision (permils)	>0.025	>0.1	0.03	0.08
Drift (permils)	>0.2	>0.8	0.13	0.33

- For water with high TDS, it takes 5-6 injections to get rid of 99% of the memory in  $\delta^{18}\text{O}$  and 98% memory in  $\delta\text{D}$ ;
- Although the Precision and Drift measurement performance of the 4% seawater solution is inferior, it remains within the specifications of the instrument.

#### 3) Memory performance evaluation with the SL

##### On DI Water



For DI water, the memory performance is stable, just as in the case when no SL is implemented.

##### On Seawater



For 4% TDS seawater samples, the memory performance starts to degrade after 200 injections (~ 30 hr of continuous run), suggesting that the SL needs to be replaced or cleaned. The replacement of the SL takes minimum down-time (a few minutes) compared with the time needed to clean the vaporizer chamber (24 hours).

By weighing the SL, we determine that 80% of the salt from the seawater solution is caught by the SL bringing the salinity of the water vapor <1000 ppm.

### Conclusions

- Precision and Drift performances are preserved for seawater isotopic analysis with the Salt Liner.
- Salt Liner effectively catches salt precipitates protecting the Vaporizer chamber from rapid salt buildup.
- Because of the decline of the memory performance caused by the salt build up in the SL, we recommend replacing and cleaning SL every 24 hr.

#### References

[1] Walker et al. (2015), Oxygen isotope measurements of seawater (H218O/H216O): A comparison of cavity-ring down spectroscopy (CRDS) and isotope ratio mass spectrometry (IRMS), *Limnol. Oceanogr.*

#### Interested in learning more?

- Contact David Kim-Hak ([dkimhak@Picarro.com](mailto:dkimhak@Picarro.com))
- Visit [www.picarro.com](http://www.picarro.com)